

XDR 3.2.11 Program Change Highlights

Version 3.2.11 (May 8th, 2017)

- Added default sharpening options and measuring for all grayscale images captured by XDR from an intraoral or pan/ceph source.
- Added 3 customizable image sharpening treatments (This may cause a slight variance from sharpening in previous versions.)
- Updated Fast Image to prevent loss of unsaved or restored images should user exit Fast Image without saving.
- Added enhancements to the Perio and Caries/Endo filters.
- Added ability to the Switch Monitor feature to be able to move between more than just 2 monitors.