

XDR Program Change Highlights

Version 3.2.15 (October 21st, 2019)

- Added Automatic Calibration File download functionality.
- Improved user feedback messages with licensing and Axium API integration.
- Expanded image recovery to include TWAIN and IOC images.
- Improved image recovery process in certain Windows 7 environments.
- Improved XDR Float to incorporate configured image treatment settings
- Improved integration of 3rd party device TWAIN sources.
- Improved exam approval function for university environments.
- Improved ability export images with image treatments.
- Improved ability to export 3rd Party saved files (e.g. pdf)
- Improved accuracy of exam date information displayed with exam.
- Added additional Sci-fi button skin selection.
- Added End User License Agreement
- New XDR Administrator Console